

PRODUCT INFORMATION

SERVA Triple Color Protein Standard II Cat. No.: 39257

PRODUCT DESCRIPTION

Application SDS PAGE Western Blotting

SERVA Triple Color Protein Standard II is a three-color, ready-to-use protein marker. It is a mixture of 10 prestained proteins covering a wide range of molecular weights (10 – 180 kDa on Tris-Glycine gels). Due to the covalently bound chromophores, 10 proteins are visible during electrophoresis or Western Blotting. The chromophore also affects the apparent molecular weight of the proteins in SDS-PAGE, unstained protein standard should be used for accurate molecular weight determination. For easy identification the 25 kDa-band is green- and the 75 kDa band is red-colored. All other protein bands are blue-colored.

Application	SDS PAGE, Western blotting	
Features	 Number of bands: 10 Three-colored Broad range size estimation: 10 -180 kDa (on Tris-Glycine gels) Easy visualization during electrophoresis and evaluation of western transfer efficiency Ready-to-use: No heating to 95 °C, dilution or addition of reducing agents before use. Supplied in a buffer containing 20 mM trisphosphate pH 7.5, 	kDa -180 -135 -100 -75 -63 -48 -35
	2 % SDS, 0.2 mM dithio- threitol, 3.6 M urea, and 15 % (v/v) glycerol.	~25
Storage	Long term storage at -20 °C: Please aliquot in a small volume to avoid repeated freezing and thawing. Stable for up to 3 months at 4 °C. Stable for up to 2 weeks at 25 °C.	711 Tris-Glycine 4~20%
Usage information	Thaw the marker at room temperature or at 3 for a few minutes to dissolve precipitates. Recommended loading volume on mini gels	7 – 40 °C

Ver. 05/17

• $2-3~\mu$ l/lane for general Western transfer For thicker or larger gels, the loading volume has to be increased.

• 3 – 5 µl/lane for clear visualization during electrophoresis

(10 x 8 cm, 0.75 or 1 mm thick):