



**One Instrument. Endless Possibilities.**  
Automated sample preparation, nucleic acid  
extraction, PCR set-up and NGS library preparation

## Flexible Solutions with Open Channel Capabilities

The *Sentosa*® SX101 is a multi-purpose liquid handling platform with open channel capabilities that automates various liquid handling steps including clinical sample extraction, PCR set-up, NGS library preparation, Sanger sequencing cleanup and many more.



### *Sentosa*® SX101 Instrument Features and Benefits


- Easy-to-implement walk away automation with pre-defined pipetting protocols for e.g. sample extraction, PCR setup and NGS library preparation
- Open channel option for customized pipetting protocols providing endless possibilities for laboratory automation
- Primary tube handling capabilities for seamless workflow integration
- Supports extraction of 21 validated clinical sample types covering a broad spectrum of sample prep requirements
- Random Access: allows assignment of different assay combinations to each sample for maximum flexibility
- Intuitive graphical user interface with easy to use touchscreen monitor
- Sample ID tracking functionality to ensure sample traceability
- Bi-directional LIS connectivity for seamless integration into existing IT infrastructure
- Continuous audit trail with user specific digital signatures in compliance to CFR 21 Part 11
- CE-IVD marked for the development of diagnostic workflows

## 21 Clinical Sample Types Validated on *Sentosa*® SX101


Sample Types (all validated)	PCR Validated Assays	NGS Validated Assays
FFPE	KRAS, NRAS, BRAF	Melanoma, NSCLC, CRC, Thyroid
Cell-Free DNA*	KRAS, NRAS, BRAF	Melanoma, NSCLC, CRC, Thyroid
Whole Blood	CMV, EBV, HSV1/2, VZV, BKV, HHV6	Leukemia*
Plasma and Serum	CMV, EBV, HSV1/2, VZV, BKV, HHV6, HBV, HCV*, ZIKV	HIV*, HCV
Cerebral Spinal Fluid	CMV, EBV, HSV1/2, VZV, BKV, HHV6	-
Urine	BKV, CMV, ZIKV	-
Stool	Rotavirus, C.Diff, Norovirus, Salmonella	-
Total RNA (Blood)	BCR-ABL m/M	-
Nasopharyngeal Swab	Flu A/B & RSV, H7N9, MERS-CoV	-
Nasal Swab	MRSA/SA, Flu A/B, RSV, H7N9, MERS-CoV	-
Genital and Oral Swab	HSV1/2	-
Sputum	MTC	-
Rectal and Perianal Swab	vanA/vanB	-
Throat and Wound Swab	Strep A	-
Lesion Swab	VZV	-
Groin, Axilla and Throat Swab	MRSA/SA	-

\*in development

## PCR



## NGS



### Real-Time PCR Workflow

In the PCR set-up, the *Sentosa*® SX101 effortlessly sample extraction and proceeds to set up samples for downstream real-time PCR in various output formats (rotor disk or 96-well).

Feature	Specifications
Throughput	8, 24 or 48 samples/ run (192 reactions/day*)  72 samples/ run (360 reactions/day*)
Time to result	~ 4 hours/run
Runs/day	4 (last run with PCR o/n)
Assay combinations	1 - 4 assays/run
Assays/day	Up to 16 different tests
Open channel capability	Yes
Reduction of system complexity	Up to 66%
Reduction of time to result by combination testing	Up to 75%
Reduction of hands-on time due to automation	20 - 30%
Reduction of service and training requirements	1 service contract
Reduction of hands-on time due to automated sample reporting	90% vs manual
Primary tube handling capability (non IVD)	Yes

\*based on 8 hour shift and assay dependant  
o/n= overnight

### Next-Generation Sequencing (NGS) Workflow

In the NGS workflow the *Sentosa*® SX101 executes the sample extraction and NGS library preparation with ease, reducing technician hands-on time.

Feature	Specifications
Throughput	Virology: Up to 15 samples/run  Oncology: Up to 7 samples/run
Time to result	~ 2 days
Open channel capability	Yes
Reduction of hands-on time due to automation	80%
Reduction of service and training requirements	1 service contract
Reduction of hands-on time due to automated sample reporting	90% vs manual

### Ordering Information

Equipment / Software	Description	Item Number
<i>Sentosa</i> ® SX101 Instrument	Automated sample processing system	400089

## DIMENSIONS

<b>Width</b>	: 107 cm (42 inches)
<b>Height</b>	: 67 cm (26 inches)
<b>Depth</b>	: 61 cm (24 inches)
<b>Weight</b>	: 95 kg (209 lb)

## ELECTRICAL REQUIREMENTS



<b>Voltage</b>	: 100 V to 240 V +/-10%
<b>Fuses</b>	: T 10 AL/250 V
<b>Frequency</b>	: 50 Hz to 60 Hz +/-5%
<b>Power Consumption</b>	: 245W
<b>Overtoltage Category</b>	: II (IEC 610 10-1)
<b>Interfaces</b>	: USB 2.0, Ethernet 100Mbit/s

## ENVIRONMENTAL REQUIREMENTS

<b>Noise Output</b>	: 56 dB (A)
<b>Heat Emission</b>	: 572 kJ/h
<b>Operating Conditions</b>	: +15 °C to + 35 °C (59 °F to 95 °F) 55% to 75% relative humidity up to 2000 m NN
<b>Storage Conditions</b>	: -20 °C to + 70 °C (-4 °F to 158 °F) 10% to 80% relative humidity

## GENERAL

<b>Protection Class</b>	: 1
<b>Degree of Contamination</b>	: 2
<b>Detector</b>	: Optical Confocal Infrared Sensor



IVD: For *in-vitro* diagnostic use. Not for distribution in US.  
Inclusive of all requisite licenses.

### Vela Operations Singapore Pte Ltd

Reg No: 201106974H  
50 Science Park Road  
#05-07 The Kendall  
Singapore 117406  
T: +65 6672 6060  
T: 1800 225 8352 (Toll-Free)  
F: +65 6672 6066  
E: infoAPAC@veladx.com

### Vela Diagnostics Germany GmbH

Albert-Einstein-Ring 15  
22761 Hamburg  
Germany  
T: +49 40 890 662233  
F: +49 40 890 662299  
E: infoEMEA@veladx.com

### Vela Diagnostics Australia Pty Ltd

PO Box 4045  
North Curl Curl 2099 NSW  
Australia  
T: 1300 012 391  
F: +61 2 9423 6970  
E: infoANZ@veladx.com

For more information, kindly contact your Vela Diagnostics representative. All rights reserved. Vela is a trademark of Vela Diagnostics Holding Pte Ltd.  
Sentosa® is a registered trademark of Vela Diagnostics Holding Pte Ltd in several markets including the US and the European Union.