

HPE[™] BlueHorizon[™]



The Better Flatbed Electrophoresis System

Best in Horizontal Electrophoresis

IEF, SDS PAGE, 2D


Broad Range of Precast Gels Available

HCP Analysis by 2D WB

IQ/OQ/PQ & Workshops on Request

Replaces MultiPhor II[™] – Compatible Format

Applications in IEF



Isoelectric Focusing

Horizontal isoelectric focusing requires high voltage. For best separation, a system with high cooling efficiency is a must as provided by the HPE[™] BlueHorizon[™] resulting in unmatched resolution. Protocols for many applications available, such as agriculture (e.g. maize, wheat, sorghum, carot, potato), veterinary (e.g. pork, cattle, turkey, chicken), food (e.g. meat, cheese, pasta) and health (e.g. MS, EPO, antibodies).

Description		Cat. no.
SERVALYT™ PRECOTES™ pH 3 - 10, size 245 mm x 125 mm x 0.3 mm	5 gels	42867.02
FocusGel 3 – 10 24S, size 250 mm x 115 mm x 0.65 mm	5 gels	43335.01
SERVA IEF Gel 3 - 9 for PhastSystem™, size 50 mm x 42 mm x 0.43 mm	10 gels	43366.01




The identification of oligoclonal bands in cerebrospinal fluid (CSF) continues to be the most sensitive diagnostic test for Multiple Sclerosis (MS). With the HPE™ BlueHorizon™ and the precast IEF gels for CSF SERVA offers a reliable system delivering higest resolution and reproducibility.

Description		Cat. no.
SERVALYT™ PRECOTES™ CSF Kit, 245 mm x 125 mm x 0.3 mm	1 kit	42800.01
FocusGel 6 – 11 40S, size 250 mm x 115 mm x 0.65 mm	5 gels	43333.01
CSF Analysis Kit for Phast™ System, 50 mm x 42 mm x 0.43 mm	1 kit	43393.01

EPO Differentiation

For characterization of natural and recombinant erythropoietin (EPO) in QC (pharma) or doping control (NADA/WADA) the EPO Doping Kit allows horizontal, IEF based sample separation. After electrophoresis the separated protein bands can be blotted onto a membrane, e.g. by pressure blotting. Secure verification of EPO isoforms is then done by immuno detection.

Description		Cat. no.
EPO Doping IEF Kit 30S	1 kit	43389.01



Blank PRECOTES[™] and Blank FocusGels

Blank PRECOTES[™] and Blank FocusGels can be equilibrated with every carrier ampholyte for highest flexibility. It allows i.a. the usage of very basic ampholytes and equilibration of the acrylamide gel matrix with urea for denaturing IEF applications.

Description	Size	Cat. no.
Blank PRECOTES™ 125 mm x 125 mm x 0.3 mm	5 gels	42759.01
Blank FocusGel 24S, 250 mm x 115 mm x 0.65 mm	4 gels	43413.01



Antibody Separation

Use SERVA precast gels and the HPE[™] BlueHorizon[™] to screen the quality of recombinant antibodies. This reproducible and reliable system can further be qualified by IQ/OQ/PQ on request.

Please contact SERVA for a customized solution at tech.service@serva.de



For cha doping tion. Aft by press

Applications in SDS PAGE

2D Electrophoresis




Horizontal 2D electrophoresis combines many advantages like higher resolution with up to 25 % more spots or much lower buffer consumption compared to traditional vertical systems. Easy handling is guaranteed due to the film-backed stabilization of the gels, available in 3 different formats.

Description		Cat. no.
2D HPE™ Large Gel NF 12.5 % Kit, 255 mm x 200 mm x 0.65 mm	1 kit	43304.01
2D HPE™ Double Gel NF 12.5 % Kit, 250 mm x 110 mm x 0.65 mm	1 kit	43302.01
2D HPE™ Triple Gel NF 12.5 % Kit, 250 mm x 110 mm x 0.65 mm	1 kit	43300.01

SDS PAGE

SDS PAGE carried out in horizontal orientation is a valuable option to replace common vertical SDS PAGE. SERVA delivers precast small and large format gels up to 255 mm x 200 mm. Efficient temperature control leads to superb resolution. Low buffer consumption, precast slots and easy handling of the film-backed precast gels are other benefits of the horizontal method.

Description		Cat. no.
SDS Gel Kit 10 % 25S, 250 mm x 125 mm x 0.45 mm	1 kit	43359.01
SDS Gel Kit NF 12.5 % 25S, 250 mm x 125 mm x 0.45 mm	1 kit	43363.01
SERVA SDS Gel 8 - 25 % for PhastSystem™, size 50 mm x 42 mm x 0.43 mm	1 kit	43503.01



SDS TA PAGE

While most precast gels for SDS PAGE contain a neutral buffer system, the Excellent gels come with a tris/acetate buffer system showing comparable band patterns like the Excel[™] gels from GE Healthcare. For separation of tissue extracts, caseins and more.

Description	Size	Cat. no.
Excellent Gel Kit 7,5 %, 250 mm x 125 mm x 0.43 mm	1 kit	43422.01
Excellent Gel Kit NF 7,5 %, 250 mm x 125 mm x 0.43 mm	1 kit	43424.01
Excellent Gel Kit 12,5 %, 250 mm x 125 mm x 0.43 mm	1 kit	43421.01
Excellent Gel Kit NF 12,5 %, 250 mm x 125 mm x 0.43 mm	1 kit	43423.01
Excellent Gel Kit 15 %, 250 mm x 125 mm x 0.43 mm	1 kit	43425.01
Excellent Gel Kit NF 15%, 250 mm x 125 mm x 0.43 mm	1 kit	43426.01

Urinary Protein Analysis

The SDS Urine Gel is optimized for urine protein separation and delivers reliable and reproducible results in horizontal gel electrophoresis. This kit allows simultaneous analysis of up to 25 samples including marker lanes, the maximum loading volume is $15 \,\mu$ l.


Description		Cat. no.
SDS Urine Gel Kit 25S, 250 mm x 125 mm x 0.45 mm	1 kit	43391.01

Antibody Separation

Apply SDS PAGE to analyse your recombinant antibodies. SERVA precast gels and HPE™ BlueHorizon™ can be qualified by IQ/OQ/PQ on request.


Please contact SERVA for a customized solution at tech.service@serva.de





Other Applications

HCP Analysis by 2D Western Blotting



The development of SERVA HPE[™] 2D Gels for Western Blotting is a powerful tool to perform state-of-the-art HCP analysis by gel-based electrophoretic separation of samples. After 1D and 2D electrophoresis, the samples are transferred to a membrane and then detected with HCP-specific antibodies.

Description		Cat. no.
2D HPE™ Triple BlotGel NF 12,5 %	1 kit	43429.01
2D HPE™ Double BlotGel NF 12,5 %	1 kit	43430.01
2D HPE [™] Large BlotGel NF 12,5 %	1 kit	43432.01

Analysis of glycans by HPE™ Electrophoresis



The high-end separation of glycans by HPE[™] flatbed electrophoresis was made possible by two parallel, interlocking innovations: (1) The development of a new chemical one pot method to release glycans from proteins including fluorescent labelling of released glycans by ANTS and (2) the development of the HPE[™] GlycoGel horizontal, film backed precast gels. The large format HPE[™] GlycoGels (250 mm x 125 mm x 0,43 mm) allow a resolution of up to 13 GUs. After electrophoresis the ANTS labelled, clearly defined glycan bands can be visualized and captured by an imaging system by applying a UV 365 transilluminator in combination with a 525 nm filter.

For more information please contact SERVA at tech.service@serva.de.



Native PAGE

Blue and clear native PAGE allow the usage of horizontal large format gels, film-backed for easy handling. Standard native PAGE and 2D native PAGE are feasible to separate large membran protein complexes.

Please contact SERVA for ordering information at tech.service@serva.de

	 	- 12 - 2		100.01	 1111	
and and	 			- Pro 1 100	 and bear	
11111	 				 	an inte
	 		1.0.0.1.1.0			

Bi-Directional Electrophoresis & 4-Electrode Lid

By using an electrode lid with three electrodes, i.e. the anode in the middle and both cathodes on the outside, the number of samples to be analysed can be doubled. This modified lid can be used whenever short runs are applicable for a reliable analysis.

A 4-electrode lid is available on request, please contact SERVA at tech.service@serva.de.

Description		Cat. no.
3-Electrode Lid for Bi-Directional Electrophoresis	1 piece	HPE-3EL



IPG Tray with Electrode Lid

Turn your HPE[™] BlueHorizon[™] into a 1st dimension IEF device for 2D electrophoresis! The HPE[™] IPG tray holds up to 12 IPG strips ranging from 7 cm up to 24 cm, each. The flexible design allows the electrode spacing to be adapted to the IPG strips used. This makes the system flexible and versatile. The high current required for IPG strip focusing is supplied by the robust 6000 V high-voltage power supply.

Description		Cat. no.
HPE™ electrode lid with IPG tray	1 kit	HPE-ELIPG
BluePower™ 6000 IPG Power Supply	1 piece	BP-6000-IPG

HPE[™] BlueHorizon[™] – optimized perfomance in flatbed gel electrophoresis

The HPE[™] BlueHorizon[™] is a flatbed system for horizontal electrophoresis using precast gels, self-cast gels and gel strips. Main applications are isoelectric focusing (IEF - like CSF, seed and food analysis, EPO, separation of antibodies and other recombinant proteins), SDS and native PAGE, 2D electrophoresis as well as separation of nucleic acids. The HPE[™] BlueHorizon[™] flatbed electrophoresis system provides unmatched resolution, reproducibility and sensitivity – the first true "High Performance Electrophoresis (HPE)" system

SERVA



Safety lid with unbreakable, robust platinum electrodes, flexible positioning

Stackable: Up to 4 units for simultaneous operation

Low buffer consumption, gels are easy to load

Drawer with safety lock, convenient to handle

SERVA

Ceramic cooling plate, easy alignment of gels

Compatible with all types of flatbed electrophoretic applications Highest resolution achieved using premium quality ceramic cooling plate Thin polyacrylamide gels, format up to 260 mm x 205 mm Broad range of precast gels for versatile applications available Complete range of reagents and consumables for self-casting of flatbed gels available IQ/OQ/PQ qualification on request Technical support and workshops Replacement of MultiPhor II[™]: compatible format

HPE[™] BlueHorizon[™] – a modular system that grows with your needs

The HPE[™] BlueHorizon[™] is a very flexible system that could easily be extended to a two-, three- or fourdeck system, still under the control of one power supply and one chiller. Stable set-up is guaranteed by stabilizing feets and clamps connecting the single units. For bi-directional electrophoresis a unit equipped with a three-electrode lid is available.



Specifications

Voltage (max.)	3000 V	
Current (max.)	25 mA	
Maximum gel size	260 x 205 mm	
Electrode distances at	270 mm, 195 mm and 115 mm	
Temperature operating range	4 °C to 30 °C	
Dimensions (WxHxD)	450 x 120 x 500 mm	
Weight	6 kg	

Ordering information

Product	Size	Cat. no.
HPE™ BlueHorizon™	1 unit	HPE-BH
BluePower™ HPE™ 3000 Power Supply	1 unit	BP-3000-HPE
HPE™ Cooling Unit (Chiller)	1 unit	HPE-CU1
HPE [™] BlueHorizon [™] PS (HPE [™] BlueHorizon [™] + BluePower [™] HPE [™] 3000 Power Supply)	1 system	HPE-BHP
HPE™ BlueHorizon™ C (HPE™ BlueHorizon™ + HPE™ Cooling Unit)	1 system	HPE-BHC
HPE [™] BlueHorizon [™] System (HPE [™] BlueHorizon [™] + BluePower [™] HPE [™] 3000 Power Supply + HPE [™] Cooling Unit)	1 system	HPE-BHSYS
HPE™ BlueHorizon™ Double Deck	1 unit	HPE-BHD
HPE™ BlueHorizon™ Triple Deck	1 unit	HPE-BHT
HPE™ BlueHorizon™ Quadra Deck	1 unit	HPE-BHQ
Stabilizing Feet, for HPE-BH	1 pair of 2	HPE-SF
Stabilizing Clamps, for HPE-BH	1 pair of 2	HPE-SC
HPE™BlueHorizon™ Bidirectional	1 unit	HPE-BH3E
3-Electrode Lid for Bi-Directional Electrophoresis	1 unit	HPE-3EL

SERVA Electrophoresis GmbH

Carl-Benz-Str. 7 69115 Heidelberg / Germany Fon: +49 6221 13840-0 Fax: +49 6221 13840-10 E-mail: info@serva.de · www.serva.de



Version 22/03